

IVD Core Materials & Solutions —

☐ ISO 13485



(b) 150 9001

## Antibody Antigen Bulk reagent for Detection



Sexually transmitted diseases (STDs) continue to pose a significant global public health challenge. Cnpair is committed to enhancing blood safety and supporting effective disease control.







Cnpair's products are designed for both 3rd and 4th generation HIV immunoassays, including reagents for anti-HIV-1/2/O and p24 antigen detection. In BBI seroconversion panel, Cnpair demonstrated 100% concordance with leading global company in CLIA.

## Comparison with Company R on BBI Panel

		Cnpair				
		Positive	Negative	Total		
Company R	Positive	30	0	30		
	Negative	0	30	30		
	Total	30	30	60		





CNPAIR BIOTECH CO., LTD.

- Mobile/ Whatsapp: +86 134 6660 2944 Haifeng Luo, PhD
- Sales: Haifeng. luo@cnpair. com
- Floor 3, Building 47, No. 488-1 Donghu North Road, Hangzhou, Zhejiang 311100, China



Hepatitis B Virus ☑ High sensitivity

☑ Mutant Detection



Cnpair provides HBsAg, HBcAb, HBeAg and HBV Pre-S1, enabling accurate identification of HBV infection. Validated by BBI panels, our reagents demonstrate superior detection sensitivity compared to other companies.

		Cnpair	R	Α	Ax	Т		
	PHM908(8)	3	3	3	3	3		
	PHM909(7)	4	3	4	3	3		
	PHM910(6)	4	4	4	-	4		
	PHM912(9)	2	2	2	2	2		
	PHM920(6)	4	4	4	4	4		
	PHM926(8)	5	5	5	5	5		
	PHM929(9)	5	5	4	5	-		
	PHM935B(12)	9	9	9	9	-		
	Total(65)	36	35	35	_	_		




Hepatitis C Virus ☑ Multiple key epitopes

☑ High accuracy



Cnpair provides chimeric and single-fragment HCV antigens, covering multiple key epitopes, enabling clear negative/positive differentiation and improving transition sample detection rate.



Syphilis

☑ Good stability ☑ Chimeric & single-fragment

Cnpair provides multiple antigen combinations, including chimeric and single fragment antigens (15/17/47KD), with reagent storage stability at 2-8°C for 15 months and a positive signal deviation of  $\leq \pm 20\%$ .

